



Scientific Solutions For Fitness™

Product Specifications

SCIFIT PRO2® All Body Exerciser

Features	<ul style="list-style-type: none">• Upper and lower crank systems individually removable, adjustable, and reversible.• Vertical seat adjustment with gas assisted lift and swivel.• Bi-directional movement and resistance.• Fit-Key™ technology, optional.• Singular lever removable seat for wheelchair access, ADA compliant.• Adjustable console angle.• Heart Rate compatible with most wireless technology.• Optional wheelchair platform with wheelchair tiedowns.• Optional high and low support boots.• Optional assist gloves and horizontal grips.
Performance	<ul style="list-style-type: none">• Resistance system-3 phase combination generator and eddy current brake with resistance output of 8 to 2000 watts.• Computer adjustable in .1 increments with 200 levels of resistance
Programs	<ul style="list-style-type: none">• Manual - One Button Quick Start.• Heart Rate - Machine adjusts workload to maintain target heart rate.• Constant Work - With wattage as a baseline, machine maintains constant workload independent of pedal speed.• Iso-Strength - for strength training.• Random - Infinite number of program profiles.• Hill Profiles - 7 hill programs with over 20 levels of intensity.
Feedback	Heart rate, time, rpm, calories, distance, level, watts, METS
Dimensions	<ul style="list-style-type: none">• Height 62"• Length 61"• Width 30"• Weight 223 lbs.• User capacity 450 lbs.
Color	Frame - Silver powder coated Covers - Metallic Gray with UV Protection
Power	Self-generating with auto recharge battery back up. Low voltage AC adaptor included.
Warranty	3 years parts, 1 year labor
Order Code	PRO230

Made in USA