

SCIFIT

PRO2[®] Total Body Exerciser

Scientific Solutions

- Bi-directional movement and resistance.
- Resistance system - 3 phase combination generator and eddy current brake with the largest wattage range in the industry, starting at 6 watts.
- 1-20 levels adjustable in .1 increments providing 200 levels of resistance.
- Oversized, swivel seat moves fore and aft, up and down and swivels in a 360 degree radius stopping at every 90 degree angle. Equipped with therapist assist levers and pedals, leaving hands free for the user
- Removable seat for wheelchair access. Accessories available.
- True adjustable step-through design, allowing all users full access and ease of entry.
- Adjustable upper and lower cranks for greater range of movement, variety and custom fit.
- All in one exerciser: Upper body, lower body or total body.
- Knee to elbow movement provides core exercise for total body wellness.
- Upper and lower cranks are dependant allowing individuals to use passive assistance.
- Heart Rate compatible with most wireless technology.

Intelli-Fit[®] Features and Programs

- 7" LCD display with user friendly interface using tactile keypad, color coded with audible beep.
- Programs: One button Quick Start, Manual, Hill Profiles, Random Profiles, Heart Rate, Iso-Strength, Constant Work, Heart-Fit Test, Heart-Fit Training (4 levels with 3 programs at each level), Power-Fit Test, Power-Fit Training (4 levels with 3 programs at each level), Stress Test and Fit-Quik[™].
- Administrative mode: includes 6 user defined languages, metric or US units of measure and other custom settings.
- Built-in twin Fans: 3 levels of intensity at low, medium and high.
- Power Burst - 8 second burst of Iso-Strength. Available on all programs except: Heart-Fit Test/Training, Power-Fit Test/Training, Heart Rate, Iso-Strength and Stress Test.
- Automatic cool down after each exercise is completed.
- Download results to USB flashdrive. Simply plug your USB key in the Fit-Key slot and your work-out will be copied in a text file for saving and printing.
- Display Feedback: Heart rate, time, RPM, calories, distance/miles, level, watts and MET's with two change view screens
- CSAFE port and 8-volt DC power port for TV or other controllers.
- Adjustable console angle to accommodate users at all positions.
- USB charger for phones, mp3 players and other electronic devices.



www.SCIFIT.com
800.278.3933

Optional Accessories



Assist Gloves



USB Fit-Keys



Wheelchair Platform
with tethers



Low Support Boots



Heart Rate Chest Strap



Wheelchair Ramp

fit-Key[®]
Compatible



Straight Hand Grips



High Support Boots

Entertainment options
also available.

Technical Specifications

- Overall Dimensions: 61" L x 30" W x 62" H
- Weight of unit: 223 lbs pre packaged
- Shipping Dimensions and weight:
68"x32"x64" 268 lbs (Skid only)
68 5/8" X 32 1/2" X 67 9/16" 303 lbs Boxed
- Maximum user weight: 450 lbs
- Self Generating with auto recharge battery backup. Low voltage AC adapter included.

Color

Frame: Silver Powder Coated
Covers: Metallic Gray with UV Protection

Warranty

3 year parts, 1 year labor. Exceptions: 1 Year warranty on pedals and seats. Wear items (rubber hand grips and pedal straps) are excluded from warranty.

Order Code

PRO230-INT